


HITACHI



Wide Open

AIRISTM II

MAGNETIC RESONANCE IMAGING SYSTEM

AIRIS™ II

*AIRIS II offers you
a premium mix of form, functionality
and performance
that are patient- and user-friendly
beyond comparison.
It's BECAUSE WE CARE.*






HITACHI

Open Technology

Ingenious open gantry design has been realized by most advanced technologies. design has been realized by most advanced technologies.




Flat gradient magnetic field coil



Flat transmitter coil



Saddle-type transmitter coil



Control panel

Open magnet circuit with two supporting columns

The gantry structure supporting the upper and lower magnets has been improved by adopting two columns in the back instead of four columns, allowing the system to provide wide open spaces of 210° in the front and 70° in the back. Needless to say, a high homogeneity of magnetic field is achieved by the optimized shimming technology. This most advanced gantry of AIRIS II has been realized only by Hitachi's unique technologies.

Flat gradient magnetic field coil


The gradient magnetic field coils are required to produce magnetic field gradients in X, Y and Z axes in the static magnetic field space. The gradient field coils used in AIRIS II are of flat type that has been developed by applying computer simulation to form a gradient magnetic field accurately with a high degree of efficiency. Because there are no resonating parts in its structure, the gradient field coils are incredibly quiet and patient friendly.

Flat transmitter coil

Development of flat-type transmitter coil is essential to realize a full-open gantry. It is necessary in the vertical field MRI system to produce horizontal high frequency magnetic field as transmitting pulse. The saddle type transmitter coil has been previously used, with which return lines exist between the upper and lower magnets, preventing the gantry to be structured as open type. Hitachi has solved this problem by developing the flat-type transmitter coil consisting of completely separated top and bottom parts. Furthermore, adoption of QD transmitter system has contributed to improvement in high uniformity and high efficiency.

Control panel

Scan operation can be started from any of the three control panels provided at the right, left and rear of the gantry, thus allowing quick examination as well as keeping communication with the patient and convenient contrast-enhanced imaging with a proper timing.



Evolutionary development toward Open Gantry

Hardware Technology

Image quality of AIRIS-II is supported by powerful hardware.

High-power gradient field power supply

Performance of the gradient magnetic field is the important basic factor to determine freedom of pulse sequence, scan time and image resolution. The high power gradient field power supply incorporated in AIRIS-II features top-class gradient field strength of 15mT/m and a slew rate of 30T/m/s. The high power gradient field strength allows thin slice imaging and high resolution imaging, while a high slew rate improves time-efficiency of pulse sequence and contributes to fast imaging.

High-output RF power amplifier

AIRIS-II incorporates RF power amplifier mounted in the rack capable of 5kW total output with very little heat dissipation. Such high power amplifier is essential to realize high performance sequences such as MTC.





Magnetic circuit with free eddy current

If the polepiece placed at the back side of the gradient field coil is made of conductive material, eddy currents are induced in the direction counteracting the gradient field, resulting in deterioration of response and degrading accuracy. AIRIS-II uses the polepiece made of nonconductive material developed to solve this problem basically.


Most powerful magnetic material

The magnet material NEOMAX used in the magnetic circuit features the most powerful coercive force. This magnet material has realized the compact, unique gantry of AIRIS-II. This permanent magnet does not require any energy at all to produce the magnetic field. Because the gantry incorporates the closed magnetic circuit, its fringe field is extremely small. The highest cost performance and installation economy are the incomparable merits of AIRIS-II. Demagnetization of the magnet material NEOMAX with time is negligibly little and its stability is outstanding.


Note: NEOMAX is the brand name registered by Sumitomo Special Metals Company, Ltd., and the material composed are Nd-B-Fe.




Deterioration of Magnet with Time



Digital RF circuit configuration



Multi-CPU Configuration



Workstation console


The operator console, a decisive factor in operability, incorporates a high performance, most advanced workstation. GUI (Graphical User Interface) operation with a mouse on a 21" color display allows users to operate the system comfortably.

Full-digital RF circuit

The performance of the RF circuit to process reception signals has a significant effect on image quality. AIRISII is equipped with a newly developed full-digital RF circuit to achieve excellent accuracy and reliability. The full-digital RF circuit incorporating a high speed AD converter digitizes reception signals in the first stage, and analog detection and filtering process which have been suffered from frequent adjustments and change with time can now be handled through data processing by the high speed digital circuits. Thus, filter characteristics can be controlled with a high accuracy.




Multi-CPU system

The 64bit multi-CPU console is equipped with the real time control pulse sequencer, dedicated high speed processing chips and image processing unit with built-in large capacity memory (512MB), thereby realizing complete multi-task processing environment to allow simultaneous operations of image acquisition, MIP processing and filming.



Advanced Technology

Development of pulse sequences without compromise AIRISII continues to provide most advanced imaging techniques and gives you answers to various clinical applications.



MR Angiography

① Multi slab

Multi slab is a technique to minimize reduction in blood flow signal in the peripheral region of the slab. In 3D-TOF sequence, the next excitation is given while blood is flowing in a slab, thereby the blood flow at the flow-out region receive excitation more than one time, resulting in reduction in signal component of the blood flow.

To avoid this, multi slabs technique in which each slab is thin and more than one slab are arrayed is used. In this technique, blood flowing in each slab in unsaturated state is excited, thus reduction in signal intensity can be minimized. With AIRISII, each slab is overlapped more than 50%, and the overlapped slices are used for image reconstruction. Furthermore, considering profile of slab, correction is made so that brightness of each slice can be flat to minimize difference between each slab. As the result, images featuring excellent continuity of blood flow and useful for diagnosis can be obtained.

Fast spin echo

② In FSE sequence image is reconstructed from measured data of more than one echo, by using multi echo, that are acquired by giving phase encode different between each echo signal. Frequency of repeatedly applied 180° pulse is called Echo Factor (EF), and scan time can be shortened by $1/EF$. A period from a 180° pulse to the next 180° pulse is called Inter Echo Time (IET). In AIRISII, these EF and IET can be changed in 1 steps according to the imaging region. Moreover, the echo positioning at the center of k-space that determines image contrast can be changed by Echo Shift (ES).

③ In FIR sequence, measurement is made in the same sequence as FSE after giving the first 180° pulse. The period from 180° pulse to 90° pulse is TI time and changing TI results in providing various image contrast. After 180° pulse excitation, longitudinal magnetization is gradually recovered, but there may be a case that the signal intensity becomes "0" in the recovery course. Even if 90° pulse is given after that, any signal cannot be obtained because no longitudinal magnetization exists. The techniques using this principle is STIR sequence for fat suppression and FLAIR sequence for CSF suppression.

④ ADA

ADA is the technique to reconstruct an image by measuring data for about a half number of phase encodes in FSE sequence and using half scan mode. For example, with EF128 k-space can be covered with one shot (one TR). In AIRISII, EF can be set up to 256, and using ADA together accelerating effect of scan time by 480 times-speed at maximum can be obtained. Application of this sequence can be endlessly widened to the clinical fields where heavily T2-weighted image is needed such as MRCP, Myelography, Urography, Ductography, etc. In these imaging applications, not only 2D measurement but also 3D measurement are usable, consequently images viewable

Neck MR angiography (2DTOF)
 Gated imaging technique is used. 2DTOF sequence causes deviation in positional information on blood flow due to heart beat, but smooth image of blood flow can be obtained by applying gated imaging technique.



Neck MR angiography (2DPC)
 Blood flow speed can be measured by gated imaging technique. These images can also be obtained simultaneously, thereby blood flow direction can be easily grasped.



Spinal stricture 68-year-old Male
 These are data of joint motion study acquired by utilizing the AIRIS-IT's feature "Open" structure. By flexing the patient forward and backward with the joint motion device available as option, the compressed status of the spine can be observed in detail. (TR/TE = 250/25)



Syringomyelia and cervical distortion 37-year-old Male
 Short DFT of 20ms is used. Because of short data sampling time, image with less flow artifact can be obtained, thereby even a small lesion cannot be overlooked. (T2-SAG/COR, TR/TE = 4000/120, RF = 11)



Cervical herniated disk
 40-year-old Female
 These are images with improved resolution by using FOV 180mm for SAG image and FOV 150mm for TRS image. A high SN ratio can be obtained with the neck solenoid coil unique to AIRISH and AMI & DWS function. Detailed information on the region compressed by herniated disk can be obtained.

Breast coil
 The matrix used either side image. The 3mm sequence in direction to observe the



Left shoulder habitual dislocation (Extinction of posterior articular lip) 20-year old
 Under AIRISH with a wide tabletop, shoulder joint can be easily positioned at the center of the magnetic field. AMI & DWS functions provide a high SN ratio with TE of in phase. (TR/TE/FA = 5000/0/0, Slice = 4.0mm)




Left frozen shoulder + periarthritis shoulder
 48-year-old Male
 The solenoid coil for joint can be flexibly set by being adjusted to the patient body. Both T1-weighted image (left) and T2-weighted image (right) can be acquired with a high SN ratio. (T1: TR/TE = 400/13, Slice = 4.0mm)



Hematoma in right femur
 T2-weighted image with flip with 15mm and 4mm slice image but is supported by clearly showing the internal (T2: TR/TE = 4000/100, RF

37-year-old Female
 Grapho-coil (option) can be
 for one-side imaging or both

co data acquired with RSSG
 can be displayed in any
 applying ADIC process to
 display.



Suspicious rerupture after meniscectomy of right knee.

33-year-old Male
 T1-weighted image with a high SN ratio can be obtained by setting very short TR to shorten scan
 time and with AMI & RWS function.
 (T1: TR/TE = 350/25, Slice = 4.0mm)



SN ratio can be obtained
 thickness in FSE sequence,
 using less number of EF, thus
 rupture of the lesion.



**Fractured media condyle &
 posterior condyle of tibia.**

37-year-old Female
 High resolution image data of thin
 slice and 512 matrix acquired by
 fully utilizing AMI, RWS and DS.
 The region necessary for diagnosis
 can be fully covered by DS.
 (T1: TR/TE = 450/20,
 Slice = 3.0mm, 16 slices)



Tumor

19-year-old Male
 With very small FOV of 150mm and
 512 matrix, sharp image of
 3mm slice is acquired.
 By selecting TE near Out of phase,
 the lesion is clearly demonstrated.
 (T2: TR/TE = 500/11.5)



Osteoncus in tibia

82-year-old Male
 Using the high sensitivity
 knee joint QD coil,
 high resolution,
 high quality image with 512
 matrix and 3mm slice
 can be acquired.
 With HET25ms and
 FOV200mm, a sufficient
 SN ratio is assured.
 (T2: TR/TE = 4000/100, EF=8)



Multiple array coil for head *



Head & neck coil *



TMJ coil *



Large joint coil *



Shoulder QD coil *



Wrist QD coil *



Breast coil *



Small joint coil *



Coils for All Regions

AIRISII offers a variety of coils with coil designs to just fit every anatomical regions.



Head QD coil



Joint coil



Flexible QD body coil
(Large-Medium) (Small#)



Knee QD coil #



Head QD coil



Head MR angiography (3DTOF)

Both multi-slice and AMI are used.

A wide region necessary for diagnosis from primary artery to peripheral arteries is imaged with a high-intensity signal. The boundaries between slices are also smooth, thus providing images for easy diagnosis.




Hexank (6) body coil (Large Myelography Special)



Myelography (3D)

High-resolution images of 3D are acquired with high EF45 and L86. By applying MP, an image quality not inferior to a conventional image without radiation can be obtained. Furthermore, by applying MIP, the region of interest with thickness can be displayed.



Multiple cerebral infarction 79-year-old Male

Signal from CSF is suppressed by FLAIR (left) sequence and small lesions of cerebral infarction near the brain surface is clearly detected. In T1-weighted (center) and T2-weighted image (right), they can hardly be identified because signal intensities in lesions and CSF are at the same level.



Suspicious vestibular schwannoma 43-year-old Female

With dual-contrast method of FSE sequence, proton density image and T2-weighted image can be simultaneously acquired. The first three echo signals are used for proton density image and the following three-echo signals are used for T2-weighted image. (TR/TE/TE2 = 4000/107/120)



Cerebral infarction sequelae

T1-weighted image (left) acquired by using both AMI and BWS. Due to short TR and TE, contrast between white matter and gray matter is excellent, and the verge of the lesions is clearly depicted.

(TR/TE = 150/15, 12 slices)



MRCP (3DFSE)

Pancreatic duct is clearly demonstrated. Measurement by ADMCP, a respiratory effect is great. With AIRIS technique can be used even in 3DFSE.



ST) areas slice, IR, image cal, reas, the vel.



Lumbar herniated disk - 27-year-old Male
Thanks to the high sensitivity flexible QD body coil and AMF & BWS function, sharp image with a small FOV can be obtained even for TRS lumbar image that requires a high SN ratio.
(T1-TRS: FOV = 200mm, TR/TE = 500/18)

5th lumbar vertebra spondylolysis and spondylolisthesis

19-year-old Male
For the region where motion artifact is little, DT 30ms is used to obtain an image with a high SN ratio. TRS image for which high resolution is required is acquired with FOV 200mm.
(T2-SAG: TR/TE = 4000/20, EF = 8)



ven for healthy ult because gated imaging quence.



Chronic pancreatitis


67-year-old Male
ADA image with IDENSE acquired by one shot. By reconstructing multiple images with AMP, the structure of pancreatic duct featuring pancreatitis can be observed.
(TR/TE = 7000/500, EF = 1.6, Slice = 3mm)



Common bile duct ectasis
69-year-old Female
Heavily hold MRCP image acquired by 2 shots using ADA. The images showing the cholangiectasis as well as pancreatic duct can be acquired in a short time with little strain to the patient. (TR/TE = 7000/500, EF = 6.4, Scan time = 14sec.)



Hysteronoma and dermoid cyst - 41-year-old Female
STIR is an effective sequence to identify the lesion as fat or others as shown in this case. Thanks to the excellent field homogeneity provided by AIRIS-II, fat signal can be suppressed evenly given for a large field of view.
(STIR: TR/TE = 3000/800/25, EF = 10)



Stone in common bile duct and occlusive cholangiopathy
50-year-old Female
Images acquired by one shot with ADA of IDENSE sequence. Stone is shown with low signal and distended bile duct can be checked in 3mm steps. (TR/TE = 7000/800, EF = 1.6, Slice = 3mm)




Liver cirrhosis
63-year-old Female
By doubling NSA using half scan technique to improve averaging effect and shortening data sampling time by setting higher BW, T1-weighted image with less artifact can be obtained.
(T1: TR/TE = 6000/16, 15 slices)

Endometrial




cyst or cyst-hyperplasia and metrorrhagia including moderate hemorrhagic component.
52-year-old Female
AIRIS-II with 0.1T magnetic field strength featuring excellent performance against artifact. Fat-suppressing nerve depressor or the like is not required for physiologically inevitable venitriculation.
(T1-SAG: TR/TE = 4000/4, Slice = 7.5mm, 14 slices)

Coil Technologies

Receiving coil is the key component of MRI system. Patient-friendly, Easy setting and High sensitivity technologies are the essential factors.



Reception signal intensity



Merits of vertical field system

AIRIS II, having two magnets at the top and bottom to generate a vertical magnetic field, can detect reception signals in the longitudinal direction. By placing a solenoid coil around the region to be imaged, sensitive regions can be effectively used. On the other hand, the conventional horizontal magnetic field systems cannot detect signals in the longitudinal direction. Therefore, saddle type or similar type coils are used as receiving coil. As shown below the saddle type coil has a structure to sandwich the imaging region with two loops, wasting a half of the sensitive zone. Difference in sensitivity between the two types has been experimentally proved to be more than 40%, that is equivalent to two times difference in scan time.

Original coil designs

AIRIS II features a rich lineup of coils. Each of them provides patient-friendliness, easy-operation and advanced technology for high image quality. The neck coil is of split type for easy setting and open type to eliminate occlusive feeling. A flexible coil fitting to the patient body can be selected. Flexible coils are designed to softly fit to any patient with a high sensitivity. Delicate considerations are given to details of easiness in operation. For example, an arm notch is provided for imaging thoracic region.

High sensitivity QD coil system




QD coil system is a technology to improve sensitivity by combining two coil systems. Detection of signals by simultaneously using two coils which have a same sensitivity region can produce two times the signals. In this case, if there is no interference between the two coils and they are independent each other, total noise level is not doubled but becomes $\sqrt{2}$. As a result, signal-to-noise ratio can be improved by about 40% higher than solenoid coil.

Multiple array coil system*

Featuring both high sensitivity and wide field of view. It is the most advanced coil system made available by high precision digital RF system. Receiving coils independent each other up to 4 channels can be simultaneously used. It is the revolution in technology that has realized both high sensitivity and wide field of view.

Principle

In general, high sensitivity and wide field of view are not satisfied at the same time in receiving coil. Therefore, a surface coil of a small diameter having a high sensitivity regionally is placed closely on the imaging region. The coil system in which several pieces of this regional coil are arrayed and used simultaneously to realize both a high sensitivity and wide field of view is the multiple array coil system. In this coil system, interference between coils is a problem. If two receiving coils are used simultaneously, they are coupled together to become equivalent to a single large coil. Such coil cannot provide a high sensitivity. This problem has been solved by the multiple array coil newly developed by Hitachi's unique technology to cancel such interference, thereby succeeded in simultaneously



Spin Echo and Gradient Echo

● AMI & BWS

AMI is the function to set asymmetric ratio of echo signals to be acquired against the peak of echo signals. Setting AMI allows improvement in contrast resolution resulting from shortened TE, increase in number of slices and improvement in S/N ratio by reduction in BW.

Band Width (BW) means frequency bandwidth of reception signal to be used for imaging. Because noise components are evenly distributed all over the frequency bandwidth, reducing BW results in reduction in noise amount to be measured. That is, in case of obtaining image with the same condition, reducing BW can improve S/N ratio of the image. In AIRISII, BW can be set to any value thanks to BWS function. Automatic setting function of AMI and BWS is also available, thereby the highest S/N ratio with any imaging parameters can be assured.

⑤ Half echo

Half echo is the function to prevent degradation in spatial resolution caused by AMI. Setting a high AMI ratio causes less number of data sampling in the high frequency region, resulting in degradation in spatial resolution.

For this reason, applying Half echo function makes image reconstruction by estimating the data for the part where data are not acquired, thereby the image resolution can be maintained without degradation.

⑥ Dual slice



Dual slice is the technique to acquire double number of slices with the same TR. The first RF pulse excites both the areas for 2 slices toward positive side, and the second RF pulse excites the first slice toward positive side and the second slice toward negative side. Adding these two kinds of data, the second slice image can be obtained by subtracting the first slice.

By using this dual slice technique together with AMI and BWS, the region necessary for diagnosis can be covered even with thin slice.

Imaging sequence	PD	T1	T2	T2*	Typical applications
Spin echo	●	●	●		General imaging
Inversion recovery		●	●		STIR image, FLAIR image
Gradient echo		●		●	Dynamic study
SARGE				●	Fluoroscopy
RS-SARGE		●			Cine imaging
TR-SARGE				●	3D imaging
Fred SE	●	●	●		Max. 480 times speed, Dual contrast, MRCP, MR myelography
Fast IR		●	●		STIR image, FLAIR image
T2* MR amplitude				●	Single/Multi slab imaging

Application Technology

Highly optimized Graphical user interface, together with the color display achieves smooth and pleasant operability.



Advanced GUI operation

GUI (Graphical User Interface) operation system, which allows the operator not only to set scan parameters but also to start scan operation and adjust display level, is employed. Since the basic operation screen is displayed with position-fixed card, the floating window does not disturb operation. Special attention is also paid to easy visualization on the 21" color display.

Easy operation

There is no limit in combination of scan parameters variable in an extremely wide range. Therefore, optimum combinations and imaging sequences can be preset according to imaging region and required image type, thereby any desired protocol set can be called and set from the clinical study library to allow easy operation.

Expert operation

In addition to the basic scan parameters such as FOV, TR, and TE, detailed parameters for AMI function, RWS function and others can also be set.

TE, which can be set in 0.1ms steps and flexible parameter setting for various high performance functions expand clinical applications of AIRIS II.

Guidance operation

If any value exceeding the specified range is entered, unique color guidance requesting attention appears, thus helps the operator to set various parameters.

Network operation

DICOM3.0[®] is incorporated in the system to support easy network operation.

Real-time operation

High speed operation on the workstation allows real-time 3D data processing.

● Real-time MIP

Unnecessary area can be deleted from 3D image data such as angiography image by delineating a free curve on projected images in three directions.

Rotating image display can also be processed in real time. Reconstruction of rotated image can be executed quickly only by setting an angle step.

● Real-time MPR

Any section image can be reconstructed in real time from 3D image data. In combination with fast 3D imaging sequence, AIRIS II shows impressive power in analysis of complicated joint and head structures.

The I-MR world that AIRIS has realized

Not only utilization of MRI system as diagnostic imaging system but also its application directly connecting to treatment. That is I-MR (Interventional MR and Intraoperative MR).

Why has I-MR become available?

Because Open type MRI systems typically represented by AIRIS series is now put in use, together with in-room monitoring function, MR fluoroscopy allowing nearly real-time imaging, open type RF receiving coil interventional MR is a clinical modality to allow the physician to lead the catheter up to the affected region using MR fluoroscopy image as guide. With the needle reached the target, various attempts have been made in combination with new treatments such as laser irradiation, injection of medicament, high frequency wave therapy, cryogenic operation, etc.

Intraoperative MR has started ahead in the field of neurosurgery treatment.

The open type MRI system is installed in the operation room and used to check immediately after craniotomy, by scanning the patient at high speed with the MRI system, if lesion such as brain tumor has been surely resected or not.



Flexible siting & Economy

MRI system helps you to explore new clinical applications.

AIRISII helps you make it.


Flexible siting

Incomparably compact design. The standard installation space is only less than 25m². The MRI scan system integrated in one unit requires no machine room.

Power consumption is only about 3kW. AIRISII features extraordinary economy. Needless to say, 5G line is within the scan room.


Economy

Compare AIRISII economy with other systems!




Specifications


Running cost



Siting cost



Cumulated expenses for 5 years



Abbreviations

- 2DFC - 2 Dimensional Phase Contrast
- 3D/SE - 3 Dimensional Fast Spin Echo
- 3DTSE - 3 Dimensional Time of Flight
- ADVA - Asymmetric Data Acquisition
- AMI - Asymmetric Measurement Imaging
- DWS - Band Width Selection
- FI-IR - Fluid Attenuated Inversion Recovery
- IET - Inter Echo Time
- MIP - Maximum Intensity Projection
- MRCP - MR Cholangio pancreatography
- MTC - Magnetization Transfer Contrast
- RSSG - RS-SARGE - RF Spoiled SAROE
- SARCE - Spontaneous Acoustic Echo Cancellation

Specifications

- **Magnetic field**
 - Static field strength 0.3 Tesla
 - ¹H resonance frequency 12.7 MHz
- **Imaging**
 - Imaging method 2-dimensional Fourier transform
3-dimensional Fourier transform
 - Imaging parameters and operations Protocol-programmed
 - Imaging coverage Whole body
 - Field of view 50 to 350 (370) mm in diameter
 - Imaging technique a) Spin echo (SE) (including Fast SE)
b) Inversion recovery (IR) (including Fast IR)
c) Gradient echo (GE)
 - Slice thickness 0.5 - 100 mm
 - Slice plane a) Transaxial
b) Sagittal
c) Coronal
d) Oblique
e) Multi-angle
 - Multi-slice 256 slices max.
 - Multi-echo Up to 4th echo
- **Measurement and reconstruction**
 - Automatic imaging function Frequency locking;
automatic tuning;
automatic gain control;
automatic shimming
 - Imaging matrix 64 to 512 X 64 to 512
 - Image reconstruction time 0.05sec. (for 256 X 256) (including display)
 - Display matrix 1024 X 1024
- **Image manipulation** Measurement of distance between 2 points;
Image profiling, ROI setting, Real-time MPR

Standard system composition

- Gantry 1
- Patient table 1
- Operator's console 1
- MRI unit 1
- ECG/peripheral pulse gating unit 1
- Filter box 1
- Standard set of accessories 1
- Receiver coils 1
 - Head QD coil (1)
 - Medium flexible QD body coil (1)
 - Large flexible QD body coil (1)
 - Joint coil (1)

Options

- Laser imager
- Knee QD coil
- Small joint coil
- Headneck coil
- Small flexible QD body coil
- DICOM 3 interface
- MR fluoroscopy
- TMJ coil
- Large joint coil
- Breast (maximum) coil
- Multiple array headneck coil
- Multiple array package
- MR angi package (TGF, PC)

Facilities

- **RF shielding facility** The scan room must be RF shielded to -60 dB in a field of 12.7±0.5 MHz in central frequency.
- **Air conditioning** Air conditioning is required to maintain the ambient room temperature and humidity within the recommended operating range. (see ENVIRONMENTAL CONDITIONS below)

Power requirements

- Single-phase 200V AC, 5kVA (for MRI system only)

Grounding facility

- 100Ω or less in grounding resistance exclusively for MRI system

Installation spaces

- **Scan room** (where the gantry and patient table are installed)
 - Standard required floor space 20 m² (≈ 4 m X 5 m)
 - Ceiling height 2.2m or more from the floor level
 - Required floor strength 3.0 ton / m² (2.5 X 3.0 m)
- **Control room** (where the MRI unit, operator's console, laser imager, air conditioner, and other components are installed)
 - Standard required space 10 m² (5 m X 2 m)

Environmental requirements

- **Ambient temperature** (a) Scan room 20 ~ 28 °C
(Temperature fluctuation: within ± 0.5 °C)

磁
遮
空
閉



• Specifications and physical appearance may be changed without prior notice for improvement of performance.
• Be sure to read Instruction Manual for correct operation of the equipment.

Hitachi Medical Corp.

Hitachi Kamakura-bashi Bldg., 1-1-14, Uchi-Kanda, Chiyoda-ku, Tokyo, Japan
International Division

Hitachi Hagoromo Bldg., 1-2-10, Uchi-Kanda, Chiyoda-ku, Tokyo, Japan